Canadian

Miscellaneous Special Publication 31 (Revised)

DFO - Library / MPO - Bibliothèque

12038814

# Fish Health Protection Regulations Manual of Compliance

626  
C314  
# 31 (REV.) D  
C.2

MISCELLANEOUS SPECIAL PUBLICATION 31 (REVISED)

QL  
626  
C314  
#31  
REV.  
D  
C.2

# Fish Health Protection Regulations Manual of Compliance

Department of Fisheries and Oceans  
Fisheries Research Directorate  
Aquaculture and Resource Development Branch  
Ottawa, Ontario K1A 0E6

Ottawa 1984

©Minister of Supply and Services Canada 1984  
Reprinted 1984

Cat. No. Fs 41-31/31-1984  
ISBN 0-662-53060-8

Available free from:

Department of Fisheries and Oceans  
Fisheries Research Directorate  
Aquaculture and Resource Development Branch  
Ottawa, Ontario K1A 0E6

Published by

Fisheries  
and Oceans

Scientific Information  
and Publications Branch

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Ottawa K1A 0E6

Correct citation for this publication:

Department of Fisheries and Oceans. 1984. Fish Health Protection  
Regulations: manual of compliance. Fish. Mar. Serv. Misc. Spec.  
Publ. 31 (Revised): 43 p.

## CONTENTS

	Page
ABSTRACT	iv
I INTRODUCTION	1
II REGULATIONS RESPECTING THE PROTECTION OF HEALTH OF FISH	3
III GUIDELINES FOR PRODUCERS	6
IV ROLE OF FISH HEALTH OFFICIALS	9
V ROLE OF LOCAL FISH HEALTH OFFICERS	11
VI SAMPLING PROCEDURES	12
VII TRANSPORTATION OF SAMPLES	16
VIII TREATMENT OF SAMPLES	17
IX PROCEDURES FOR THE DETECTION OF CERTAIN BACTERIAL FISH PATHOGENS	18
X PROCEDURES FOR THE DETECTION OF VIRUSES	24
XI PROCEDURES FOR THE DETECTION OF CERTAIN PARASITES	28
XII EGG DISINFECTION PROCEDURES	31
REFERENCES	33
APPENDIX 1 NATIONAL REGISTRY OF FISH DISEASES	35
APPENDIX 2 REGIONAL ADMINISTRATIVE AUTHORITIES	36
APPENDIX 3 QUALIFICATIONS OF FISH HEALTH OFFICIALS	39
APPENDIX 4 FISH HEALTH CERTIFICATE	40
APPENDIX 5 FISH HEALTH LABORATORY REPORT	42

## ABSTRACT

Department of Fisheries and Oceans. 1984. Fish Health Protection Regulations: manual of compliance. Fish. Mar. Serv. Misc. Spec. Publ. 31 (Revised): 43 p.

This revised manual explains the application of the Fish Health Protection Regulations under the Fisheries Act of Canada. It outlines the administrative and inspection procedures to be followed and provides step-by-step procedures for handling fish samples to test for the major bacterial, viral, and myxosporean pathogens of salmonids.

All movements of fish into Canada or between provinces require a permit. This permit will be issued only to producers whose facilities have been inspected and certified free of designated disease agents.

The sampling procedures are based on the probability of detecting an infected specimen in a lot, assuming a certain prevalence of detectable infection. Selection, transportation, and laboratory handling of the samples are described in detail.

The methods provide for the detection of the following pathogens: the kidney disease bacterium (Renibacterium salmoninarum); the redmouth bacterium (Yersinia ruckeri); the furunculosis bacterium (Aeromonas salmonicida); the protozoans causing whirling disease (Myxosoma cerebralis) and ceratomyxosis (Ceratomyxa shasta); the viruses causing viral hemorrhagic septicemia, infectious hematopoietic necrosis, and infectious pancreatic necrosis; and other pathogens considered to be notifiable.

## I INTRODUCTION

The future of Canadian fish culture and of recreational and commercial fisheries depends upon healthy fish stocks. International and interprovincial exchange of cultivated species is increasing rapidly, facilitating the dissemination of serious infectious diseases, for example the spread of whirling disease in North America. Effective programs of prevention and control are necessary to prevent spread of fish disease agents and to alleviate this serious impediment to development of these fisheries. This revised manual of compliance was prepared to explain the application of the Fish Health Protection Regulations under the Fisheries Act of Canada. It presents guidelines for producers, defines the roles of Fish Health Officials and Local Fish Health Officers, and outlines the sampling, handling, and diagnostic procedures that constitute inspections leading to certification.

The Fish Health Protection Regulations as published in the Canada Gazette appear in Section II. They apply to all fish species belonging to the family Salmonidae of the genera listed in Schedule I. The Regulations are designed to prevent the spread of infectious diseases through inspection of production facilities' fish stocks, and to control the movement of infected fish. They apply to live and dead cultured fish, eggs (including any fertilized or unfertilized sex products) of cultured and wild fish and products of dead cultured fish (including cold smoked products, but not heat processed canned products) destined to move into Canada or across provincial boundaries within Canada. In the event of violation of these Regulations, seizure and other powers of the Fisheries Act apply.

Compliance with these Regulations requires that four consecutive satisfactory inspections over a period of not less than 18 months must be conducted at a facility prior to certification. Inspections must be conducted at intervals of not less than 90 days and not more than 270 days during the pre-certification period. A new facility using an isolated water supply free of all species of fish, starting with stocks from a certified facility, can gain certified status after only one satisfactory inspection. For a production facility to retain certified status, consecutive satisfactory inspections must be conducted twice yearly at approximately 6-month intervals (not less than 90 days nor more than 270 days). Inspections should be made during the spring and fall periods, with the final inspection having been made not more than 270 days prior to date of entry of a shipment of fish into a province. Procedures outlined in this manual are to be followed.

The concept of "zoning" may be introduced on the basis of detailed information relating to the absence of disease agents in specific areas. Recognition of specific pathogen free "zones" could lead to modification of inspection requirements in production facilities. The zoning concept will only be applied to facilities exporting dead cultured fish.

If further clarification on the interpretation of these regulations is required, contact the National Registry of Fish Diseases, Department of Fisheries and Oceans, Ottawa, Ontario K1A 0E6, telephone (613) 990-0276 or telex 0534228.

## II REGULATIONS RESPECTING THE PROTECTION OF HEALTH OF FISH

### Short Title

1. These Regulations may be cited as the Fish Health Protection Regulations.

### Interpretation

2. In these Regulations,  
"approved" means approved by the Minister; (approuvé)  
"certificate" means a certificate issued in accordance with section 6; (certificat)  
"cultured fish" means a fish listed in Schedule I that is propagated by man in a fish culture facility and includes the eggs of such fish; (poisson d'élevage)  
"fish health official" means a person approved to inspect fish and fish sources for the purposes of these Regulations; (inspecteur sanitaire des poissons)  
"import" means to bring into any province of Canada from any other country or any other province of Canada; (importer)  
"import permit" means a permit issued pursuant to section 4; (licence d'importateur)  
"local fish health officer" means a person approved as a local fish health officer in charge of the administration and enforcement of these Regulations; (agent local de protection de la santé du poisson)  
"Minister" means the Minister of Fisheries for Canada; (Ministre)  
"wild fish" means a fish listed in Schedule I, other than a fish propagated by man in a fish culture facility. (poisson sauvage)

### Prohibition

3. No person shall import cultured fish or eggs of wild fish without an import permit.

### Permits

4. Subject to section 5, a local fish health officer for a province may issue an import permit to a person who applies therefor authorizing him to import cultured fish or eggs of wild fish into that province.
5. No import permit shall be issued unless
  - (a) the person who applies for the permit has obtained a certificate; and
  - (b) the local fish health officer is satisfied that any disease or disease agent listed in Schedule IV indicated on the

certificate pursuant to subsection 6(2) will not be harmful to the protection or conservation of fish in the province of importation.

### Certificates

6. (1) A certificate required by paragraph 5 (a) shall be issued by a fish health official certifying in respect of
- (a) live cultured fish that
    - (i) the fish come from a source that was inspected in an approved manner and found to be free of any disease or disease agent listed in Schedule II, and
    - (ii) no fish, other than a fish from a source referred to in subparagraph (i), has been introduced to the source of importation within the two years immediately preceding the date of certification;
  - (b) eggs of wild fish, that the eggs were taken from wild fish that were inspected in an approved manner and found to be free of any disease or disease agent listed in Schedule II; or
  - (c) dead cultured fish, that the fish come from a source that was inspected in an approved manner and found to be free of any disease or disease agent listed in Schedule III.
- (2) A certificate referred to in paragraph (1)(a), in respect of a live cultured fish, or paragraph (1)(b) in respect of eggs of wild fish, shall have indicated thereon the presence of any disease or disease agent listed in Schedule IV that has been detected in the fish during the inspection referred to in that paragraph.

### SCHEDULE I

#### FISH

All species and hybrids derived from species of fish belonging to the family Salmonidae, including the following genera:

Pacific salmon	<u>Oncorhynchus</u> spp.
Danube salmon and Taimens	<u>Hucho</u> spp.
Atlantic salmon	<u>Salmo</u> spp.
Trout	<u>Salmo</u> spp.
Char	<u>Salvelinus</u> spp.
Grayling	<u>Thymallus</u> spp.
Lenok	<u>Brachymystax</u> spp.
Inconnu	<u>Stenodus</u> spp.
Whitefish	<u>Coregonus</u> spp.
Whitefish	<u>Prosopium</u> spp.
Ayu	<u>Plecoglossus</u> spp.

## SCHEDULE II

### DISEASES OR DISEASE AGENTS FOUND IN LIVE FISH OR THEIR SOURCE

1. Any filterable replicating agent capable of causing cytopathic effects in the cell lines of fish specified by the Minister, including but not limited to
  - (a) Viral Hemorrhagic Septicemia (Egtved) (Egtved virus, VHSV)
  - (b) Infectious Hematopoietic Necrosis (IHNV)
  - (c) Infectious Pancreatic Necrosis (IPNV)
2. Whirling Disease (Myxosoma cerebralis)
3. Ceratomyxosis (Ceratomyxa shasta)
4. Furunculosis (Aeromonas salmonicida)
5. Bacterial Kidney Disease (KD bacterium)
6. Enteric Redmouth Disease (RM bacterium)

## SCHEDULE III

### DISEASES OR DISEASE AGENTS FOUND IN DEAD FISH OR THEIR SOURCE

1. Viral Hemorrhagic Septicemia (Egtved) (Egtved virus, VHSV)
2. Whirling Disease (Myxosoma cerebralis)

## SCHEDULE IV

### DISEASES OR DISEASE AGENTS FOUND IN LIVE FISH OR THEIR SOURCE

1. Myxobacterial infections
2. Motile Aeromonad Septicemia (motile Aeromonas sp.)
3. Pseudomonad Septicemia (Pseudomonas spp.)
4. Vibriosis (Vibrio spp.)

### III GUIDELINES FOR PRODUCERS

The Regulations apply to all facilities from which live and dead cultured fish, eggs of cultured and wild fish, and products of dead cultured fish (including cold smoked products but not heat processed canned products) of all species belonging to the family Salmonidae, as designated in Schedule I, which will be shipped into Canada or from one province to another. "Facilities" include all sites used for the propagation or holding of eggs and/or fish.

An import permit is required for each shipment of fish or eggs designated above, except for interprovincial shipments of dead cultured fish in Canada. These permits may be issued only to those producers providing appropriate evidence that the health status of their facility has been certified through inspection of their fish stocks by an approved Fish Health Official. Producers should become familiarized with all of the regulations of the region into which they plan to import fish since requirements may exist in addition to the Fish Health Protection Regulations.

A list of Fish Health Officials approved for the certification of facilities may be obtained from the National Registry of Fish Diseases, Department of Fisheries and Oceans; however, Canadian producers should request certification directly from the authority under which they are licenced to operate.

Within Canada, a producer wishing to obtain certification for his/her production facility is expected to provide fish at his/her own expense for inspection. Fish species, other than those listed in Schedule I, that the producer is rearing at his/her facility are also subject to sampling by the Fish Health Official. The Department of Fisheries and Oceans will assume the costs of such inspections within Canada, only as they relate to commercial operators who ship fish outside the province in which they are licenced. Outside Canada, the producer must obtain certification at his/her own expense using Fish Health Officials approved by the Minister. The Fish Health Official has flexibility in establishing time and frequency of sampling and selection of fish and must have access to records relating to introductions, losses, disease prevalence and treatments.

If a production facility which is certified receives fish or eggs from an uncertified source including species not listed in Schedule I, or if any disease specified in Schedule II in the Regulations is found in the facility, it is subject to de-certification and a 2-year re-inspection program prior to re-certification. Before commencing a re-certification program, the producer must demonstrate that he/she has undertaken adequate measures to eliminate actual or potential disease problems. If a certified facility receives fish or eggs from another certified source, the certification status is unaffected, but that producer must be prepared to provide a record of such introductions upon request. A new production facility located on an isolated water supply free of all species of fish may be certified after only one satisfactory inspection if the stock(s) in that facility was obtained from another certified facility. After this first inspection, maintenance of certification status will depend on continuing

satisfactory inspections similar to that required by all other facilities.

Before shipment, an import permit must be obtained from the Local Fish Health Officer of the province to which the shipment is going. Import permits must accompany all shipments of fish. Offices at which the Local Fish Health Officers can be reached are listed in Appendix 2.

To acquire an import permit, producers must provide a signed copy of their Fish Health Certificate (Appendix 4) attesting that no fish, other than those from certified facilities or certified wild stocks, including species not listed in Schedule I, have been introduced into their facilities during the period of certified status. When an application for an import permit is made, the producer must also declare that no certifiable disease agent has been found at the facility since the last inspection.

An import permit allows the transport of fish from the certified source facility directly to the receiving facility. The replenishing of water in a truck or other shipment container at any location other than a certified facility or an isolated water source free of all species of fish is the same as an introduction of fish from an uncertified source into a facility, and invalidates the import permit.

Although not specified in the Regulations, surface disinfection of eggs prior to shipment is strongly recommended unless otherwise advised. A suggested procedure for disinfection is given in Section XII of this manual.

Specialized producers may need clarification on some topics. Shipments of fish distributed dead from landlocked lakes and prairie lakes (potholes) that regularly undergo winter kill are a special case: they may be shipped without inspection if (a) the lakes were stocked from a certified source, (b) no overt infectious disease has occurred during the growing period, and (c) the lakes have no history of contamination with the certifiable disease agents listed in Schedule III. If all of the foregoing conditions cannot be met, the lakes must be certified as being free of the disease agents listed in Schedule III before shipments of dead fish may be made. If live fish are to be shipped, the lakes must be certified as being free of the disease agents listed in Schedule II. Certification must be based on the outcome of inspections of fish produced in the lakes conducted at, or about, the time of fish harvesting.

A Fish Health Official may wish to alter inspection procedures under certain circumstances. A prior ruling on these alterations can be obtained from the National Registry of Fish Diseases.

A fish health diagnostic laboratory may be located in a different province from that in which an inspection is being conducted, requiring the transport of fish samples across provincial boundaries.

Such samples should be handled so as to avoid any transfer of disease. All persons involved in fisheries research must conform to the requirements and intent of the Regulations and be cognizant of regional conditions.

#### IV ROLE OF FISH HEALTH OFFICIALS

The Fish Health Official must be a qualified specialist in fish disease diagnosis, have access to a laboratory equipped to undertake the diagnostic procedures outlined in this manual, and have received approval from the Government of Canada (Appendix 3). Those seeking approval are to contact the National Registry of Fish Diseases, providing a copy of their curriculum vitae, three specimens of their signature, and an outline of their laboratory capabilities. All applications will be screened by a panel of disease specialists, and applicants will be advised as to the outcome of their application. Evidence that an approved Fish Health Official is not following the intent of the Regulations, i.e. to prevent the spread of infectious disease agents by careful inspection of production sources and control of infected stock movements, by wilfully disregarding either inspection requirements and/or results of inspections, will lead to suspension or removal from the approved list.

The Fish Health Official conducts an inspection by visiting the production site, viewing all parts of the site, and carrying out the procedures outlined in Sections VI-XI of this manual. The Official should obtain information from the owner or manager of the site on identification of stocks being inspected and familiarize himself/herself with records of introductions, losses, disease prevalence and treatment of fish in the facility over the past two years, or from the date of initial certification. If all requirements of the Regulations are met and if the inspection schedules outlined earlier have been adhered to, the Fish Health Official may issue a Fish Health Certificate (Appendix 4), which must be distributed as follows:

- White copy - to owner/manager, a copy of which will accompany each application for an import permit
- Yellow copy - to National Registry of Fish Diseases
- Pink copy - to Fish Health Official for file use

If Fish Health Officials have diagnostic information which they have gathered from sampling outside of the designated spring-fall sampling periods, or documented diagnostic information from other reliable sources, this information is to be used in determining certification status along with that collected during the designated periods.

In addition, Fish Health Officials in Canada are required to complete a Fish Health Laboratory Report (Appendix 5) in triplicate, sending the yellow copy to the National Registry of Fish Diseases, the white copy to owner/manager to provide evidence of the inspection, and retaining the pink copy for the files of the Fish Health Official. Fish Health Laboratory Reports of facilities outside Canada do not have to be filed with the Registry, although the laboratory results leading to certification of a facility may be required by the purchaser.

If de-certification of a facility is necessary because of failure to meet the various requirements of the inspection, either through detection of Schedule II pathogens or through introduction of stocks from a non-certified source, then the following procedures are to be carried out by the Fish Health Official: copies of the completed Fish Health Laboratory Report and a letter outlining the reason(s) for de-certification and the steps necessary to regain certification are to be sent to the producer and the National Registry of Fish Diseases.

The minimal steps necessary to regain certification of a facility are as follows:

- a) All lots in which a Schedule II disease agent was detected must be removed from the facility.
- b) The rearing units from which these lots were removed must be treated using recognized procedures to eliminate disease agents.
- c) The facility must undergo four consecutive satisfactory inspections over a period of not less than 18 months, with intervals of not less than 90 days and not more than 270 days between inspections.

## V ROLE OF LOCAL FISH HEALTH OFFICERS

Local Fish Health Officers, located within each region of Canada, administer these Regulations for their province or area. Their responsibilities include a review of certificates and data relevant to a production source in question, a shipment of fish or eggs in particular, and the health needs of their region. Subsequent to Section 5 of the Regulations, they may issue import permits to applicants to allow passage of acceptable shipments of live fish or eggs entering Canada or moving interprovincially within Canada and dead fish entering Canada at border points. Permits must accompany each shipment. A standard import permit is not provided with this manual since provinces or regions have developed import permits in response to their individual needs.

The notifiable disease agents in Schedule IV are those not listed in Schedules II or III, but considered important by the Fish Health Official. In special cases, just cause may exist for the Local Fish Health Officer to prevent the import of these notifiable disease agents or other fish pathogens into an area.

Local Fish Health Officers are located at offices listed in Appendix 2.

## VI SAMPLING PROCEDURES

### A. CULTURED FISH

#### 1. Sampling by lot

Unless noted (see VI B 1), sampling of fish shall be according to "lot." A lot is defined as fish of the same age that have always shared the same water supply and that have originated from a discrete spawning population. In situations where this lot definition cannot be applied, the Fish Health Official will use his/her own discretion in dividing the fish into lots.

#### 2. Selecting the sample

The method of determining the actual number of fish to be sampled from any particular lot is based on obtaining a 95% probability of detecting an infected specimen in a lot with an assumed prevalence of detectable infection of 5 or 10% (Table 1). It is important to note that certain disease agents when in the carrier state are very difficult to detect. The Fish Health Official and the producer must be aware of this possibility. Statistical probabilities as stated in Table 1 may not apply in such situations. The samples must be taken under the supervision of the Fish Health Official in a manner and at a time that tends to provide the greatest opportunity of detecting any disease agent present. When the sample is being withdrawn from a lot held in a single holding unit (e.g. tank, raceway, or pond), the sample must be selected to contain as many moribund and freshly dead specimens as are available. If the previously mentioned lot of fish is being held in more than one holding unit, the total

TABLE 1. Sample size required to detect one or more infected specimens in populations (lots) with an assumed minimum prevalence of detectable infection of 5 and 10%. Calculations are based on a 95% level of confidence. For intermediate population sizes, use the sample size for the next larger population listed. (From Ossiander and Wedemeyer, 1973.)

Population size	Number of fish to be sampled when assumed prevalence of detectable infection is:	
	5%	10%
50	29	20
100	43	23
250	49	25
500	54	26
1,000	55	27
2,500	56	27
5,000	57	27
10,000	57	27
100,000	57	27
Over 100,000	60	30

number of specimens to be collected would be the same as before. However, the specimens must be taken from the individual units in numbers that reflect the proportion of the lot held in each of the units. Again, the sample from any given unit must consist of as many moribund and freshly dead fish as are available. If additional fish are required to complete the sample, apparently healthy fish may then be collected, or the Fish Health Official may elect to return within a 30-day period to complete collection of the required number of fish.

The samples should not be taken during or immediately after therapeutic treatment. Complete background information for all samples must be obtained. This includes data relating to any recent use of chemotherapeutants, the health history of the facility, and the lots from which the samples were taken.

In the event that overt disease signs are noted at the time of sampling, procedures for the detection of notifiable disease agents and other pathogens should be carried out in addition to the procedures for identification of disease agents listed in Schedule II.

Where samples of fish or their tissues are to be processed in pools, as in the virological assays, rather than individually, care must be taken to process those from apparently healthy fish separately from those of freshly dead or moribund fish.

### 3. Sample size for viruses

- a) Production fish (non-broodstock fish): the sample size to be employed is one that gives a 95% probability of detecting an infected specimen in a lot, assuming the minimum prevalence of detectable infection is 5%.
- b) Broodstock (sexually mature fish held for reproductive purposes): the same sensitivity in detecting infected individuals must be aimed at. Sampling must be conducted once yearly at spawning time. In those species for which spawning is a terminal event, tissue samples must be collected from all fish involved up to a maximum of 60 fish. For those species that are repeat spawners, 10%<sup>1</sup> of all spawners used up to a maximum of 30 fish, must be subjected to lethal sampling. The balance of the samples required to achieve the rate that provides a 95% probability of detecting an infected fish in a lot, assuming the minimum prevalence of detectable infection is 5%, is to be made up of reproductive fluids. Ovarian fluid must account for as many as possible of the reproductive product samples collected.

---

<sup>1</sup>The lethal sampling of only 10% of all the spawners does not conform to the sampling rates set out in Table 1; it is intended to conserve small but valuable populations of spawners that may live to spawn again.

4. Sample size for bacteria

- a) Production fish must be sampled for bacterial pathogens at a rate that provides a 95% probability of detecting an infected specimen in a lot assuming the minimum prevalence of detectable infection is 5%. For routine purposes, only fish averaging four cm or more in fork length need be bacteriologically sampled. Reliable bacteriological sampling of fish smaller than four cm in fork length is technically more difficult; such fish need be sampled only when unusual and unexplained mortality rates or disease signs are observed.
- b) Broodstock must be sampled for bacteria in accordance with the sample size already outlined for the lethal sampling of broodstock for virus (see VI A 3 b).

5. Sample size for parasites

- a) Production fish must be sampled for Myxosoma cerebralis and Ceratomyxa shasta at a rate that provides a 95% probability that an infected fish will be detected in a lot assuming the minimum prevalence of detectable infection is 5%.  
M. cerebralis: the fish must be at least 120 days of age for the test to be meaningful because the spores on which the diagnosis is based are slow to develop. In fish grown or held at temperatures below 12°C, spore formation may take 9-11 months (Taylor et al. 1973).  
C. shasta: routine monitoring of apparently healthy fish for the presence of spores need be performed only with fish averaging at least 120 days of age. Younger fish need be examined for C. shasta only when unusual and unexplained mortality rates or disease signs are observed.
- b) Broodstock must be sampled in accordance with the sample size already outlined for the lethal sampling of broodstock for virus (see VI A 3 b).

6. Times of sampling and frequency

- a) Production fish: sampling will be conducted at least twice yearly with sampling times and frequency depending on local conditions and the discretion of the Fish Health Official. Because the detection of M. cerebralis and certain viruses is best accomplished using fish of a certain age, spring and fall sampling periods (March-May and September-November) are recommended.
- b) Broodstock: sampling will be conducted once yearly at spawning time (see VI A 3 b).

B. WILD FISH

1. Sexually immature fish

All sexually immature fish taken from the wild must be sampled at a rate that will give a 95% probability of detecting an infected specimen in the total catch, assuming the minimum prevalence is 5%. This sampling rate will apply to the types of pathogens (viral, bacterial, and myxosporean) mentioned earlier. The Fish Health Official may use considerable discretion when determining an appropriate sample size for wild fish.

2. Sexually mature fish

If wild spawners or their fertilized eggs are to be collected, seminal and ovarian fluids must be taken from all the fish involved, up to a maximum of 60 fish. In those species for which spawning is a terminal event, tissue samples must be collected from all fish involved (lethal sampling) up to a maximum of 60 fish. For those species that are repeat spawners, 10%<sup>2</sup> of all spawners used or collected, up to a maximum of 30 fish, must be subjected to lethal sampling. These sampling rates apply to the types of pathogens (viral, bacterial and myxosporean) mentioned earlier. Again, considerable discretion must be used by the Fish Health Official when certifying wild stocks.

C. FERTILIZED EGGS AND SEX PRODUCTS

The sampling of fertilized eggs or the sex products of fish cannot be relied on to detect disease agents listed in Schedule II. The threat to fish health posed by such eggs or sex products must, therefore, be determined by examining the disease status of the parent fish.

D. NON-SALMONID SPECIES

Other species not belonging to the family Salmonidae occurring at the same facility as salmonid species are subject to sampling by the Fish Health Official for disease agents listed in Schedules II or III.

---

<sup>2</sup>See footnote 1.

## VII TRANSPORTATION OF SAMPLES

Fish samples must be handled rapidly in such a way that degenerative changes do not render diagnosis unreliable or impossible. If samples cannot be brought to the laboratory alive, they should be stored on ice or refrigerated for no longer than 48 hours.

### A. LIVE FISH

Live fish should be transported in sealed plastic bags that have been partly filled with water and charged with oxygen. The bags may then be placed with ice in insulated containers. Under these conditions anesthetics are not usually required.

### B. DEAD FISH

Fish sampled should be placed in sealed plastic bags (dead or moribund fish should be kept separately from healthy fish) that are packed in an insulated container with a layer of ice around each bag.

### C. REPRODUCTIVE FLUIDS

Seminal and ovarian fluid samples should be shipped in sealed containers on ice. Do not mix seminal and ovarian fluid samples. For pooling restrictions see X B 4 b and X C.

## VIII TREATMENT OF SAMPLES

### A. AUTOPSY PROCEDURE

#### 1. General comments

The procedure outlined below is designed to facilitate the processing of large numbers of fish for the designated pathogens. Except for very small specimens, the same fish will serve as the source of tissues for the various bacteriological, virological, and myxosporean tests. The bacteriological examination must be performed first. To maximize detection sensitivity, fish must be autopsied within 48 h of sampling and all assay procedures should be initiated within this time.

#### 2. External examination and sampling

Note all gross abnormalities such as body discoloration, body distention, exophthalmia, ulcers, blebs, inflammation, hemorrhagic areas, clubbed or necrotic gills, and eroded opercula, fins and caudal peduncles. Inoculate the appropriate media and prepare stained smears with material from these lesions.

#### 3. Internal examination and sampling

Disinfect the surface of the fish and using aseptic techniques expose the kidney. Inoculate the appropriate media and prepare the appropriate smears with kidney material.

Inspect the viscera for abnormalities. Inoculate the appropriate media and prepare smears with material from any abnormal organs. Remove tissues for virological and myxosporean examinations.

### B. DISPOSAL OF SAMPLES

The receiving laboratory should handle and dispose of samples and other items likely to be infectious in a manner that precludes the dissemination of disease agents. All materials such as fish carcasses or tissues, transport containers and water, microbial cultures, and contaminated equipment should be autoclaved, incinerated, or otherwise sterilized before being discarded.

**IX PROCEDURES FOR THE DETECTION OF CERTAIN  
BACTERIAL FISH PATHOGENS**

**A. SCOPE**

The microorganisms covered by these procedures are divided into three categories:

1. Those pathogens (listed in Schedule II) considered to be limited in geographical distribution and the absence of which must be verified. These include the following certifiable disease agents:  
Renibacterium salmoninarum (Bacterial kidney disease)  
Yersinia ruckeri (Enteric redmouth disease)  
Aeromonas salmonicida (Furunculosis)
2. Those pathogens (listed in Schedule IV) that may be ubiquitous, and if detected during a disease outbreak or systemically in the absence of clinical signs in sampled fish, are notifiable. These include:  
Myxobacteria (e.g., Flexibacter columnaris)  
Motile Aeromonas sp. (e.g., Aeromonas hydrophila)  
Pseudomonas spp. (e.g., Pseudomonas fluorescens)  
Vibrio spp. (e.g., Vibrio anguillarum)
3. Those bacterial agents not listed in Schedules II and IV, which the Fish Health Official detects and determines to be associated with significant losses and/or disease signs (e.g., Edwardsiella tarda).

**B. PROCEDURES**

The following procedures represent the minimum requirements in bacteriological testing; they must be applied to all samples taken from lots in which there is an unusually high prevalence of disease signs and/or mortalities. In samples from apparently healthy lots the procedures need only be applied to those fish averaging four cm or more in fork length (see VI A 4 a).

1. Aseptically obtain the following tissues/material and streak on the appropriate medium:
  - a) kidney tissue, preferably from sites that appear abnormal, and external/internal lesion material on Tryptic Soy Agar (TSA);
  - b) gill tissue and/or external lesion material on Shieh's medium (SH agar) or Cytophaga Agar (CA) only if gross pathology suggests a myxobacterial infection.
2. Prepare Gram-stained smears of kidney tissue and lesion material and examine a minimum of 25 fields (900-1000X magnification). The presence of small Gram-positive diplobacilli, frequently present intracellularly, is presumptive evidence for R.

# ERRATA

Department of Fisheries and Oceans. 1984. Fish Health Protection Regulations: manual of compliance.  
Fish. Mar. Serv. Misc. Spec. Publ. 31 (Revised): 43 p.

**The following figure is a replacement for the one found on page 19.**


FIG. 1. Differentiation of bacterial isolates from fish on TSA, SH or CA plates.


FIG. 1. Differentiation of bacterial isolates from fish on TSA, SH or CA plates.

salmoninarum (Sanders and Fryer 1980). (Presumptive evidence is strengthened by lack of growth on TSA.)

3. Incubate the TSA plates at 20°C for five days; examine daily for growth. Incubate the SH or CA plates at 15-20°C for five days; examine daily for growth.
4. If gross pathology suggests a myxobacterial infection, prepare moist mounts of gill tissue and/or lesion material and examine for the presence of masses of long, thin rods. Select young representative colonies (yellow dry, rhizoid or yellow, moist, spreading) from the SH or CA plates and prepare Gram-stained smears. The presence of long, thin, Gram-negative rods, capable of gliding or creeping motility, constitutes a presumptive diagnosis for myxobacteria.

5. Select young representative colonies from TSA plates. Differentiate the microorganisms on the basis of the following characteristics:

- a) If the cells are Gram-negative, rod-shaped, oxidase-positive, non-motile, ferment glucose (O.F. test) and usually produce a brown diffusing pigment, the isolate is presumptively A. salmonicida (Griffin et al. 1952). Achromogenic strains of A. salmonicida may occur (Evelyn 1971).
- b) If the cells are Gram-negative, rod-shaped, oxidase, indole and H<sub>2</sub>S-negative and produce an alkaline/acid (K/A) reaction in Triple Sugar Iron Agar (TSI), the isolate is presumptively Y. ruckeri.
- c) If the isolate differs from 5 b) by being indole and H<sub>2</sub>S-positive and produces both acid and gas in TSI, it is presumptively E. tarda (Amandi et al. 1982).

A flow chart is given for these procedures in which additional features separating motile aeromonads, Pseudomonas spp., and Vibrio spp. are indicated (Fig. 1).

Confirmatory testing of presumptively identified certifiable agents must be performed. For a possible exception, see IX C 2 g. The serological and biochemical tests described in IX C 2 are to be used for confirmatory testing.

Confirmatory tests for the notifiable agents may be performed at the discretion of the Fish Health Official. Identification should be based on whether or not the disease agent is associated with significant losses and/or gross pathology.

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Agar	1.0%
Peptone	0.5%
Yeast extract	0.05%
Magnesium sulfate	0.03%
Sodium pyruvate	0.01%
Monobasic potassium phosphate	0.01%
Dibasic potassium phosphate	0.005%
Sodium bicarbonate	0.005%
Calcium chloride	0.001%
Citric acid	0.001%
Sodium acetate	0.001%

---

<sup>3</sup>The products specified have proven satisfactory for the purposes indicated; this, however, does not imply that other products may not be equally satisfactory.

Barium chloride	0.001%
Ferrous sulfate	0.0001%

pH	7.0
----	-----

NOTE: This is a modification of Shieh's medium. Neomycin (5 µg/mL) and polymyxin B (10 units/mL) can be added to SH agar to facilitate the isolation of myxobacteria by suppressing the growth of other bacteria (Fijan 1969).

c) Cytophaga Agar (Anacker and Ordal 1959)			
Tryptone	0.05%	Beef extract	0.02%
Yeast extract	0.05%	Agar	0.9%
Sodium acetate	0.02%	pH	7.2-7.4

## 2. Test media, reagents, and methods

- a) Cytochrome oxidase test  
Using a platinum loop, transfer some bacterial growth from an actively growing plate culture to a paper strip impregnated with the appropriate chemicals. After thorough spreading, a positive test is indicated by the development of a bright blue color within a minute (MacFaddin 1980).
- b) Motility  
Examine log-phase cultures in wet preparations using Tryptic Soy Broth as the suspending medium. If the wet preparation method gives a doubtful result, check the results by stab-inoculating tubes of Motility Test Medium (Difco) or Glucose Motility Deep (GMD) medium (Walters and Plumb 1978).
- c) Differentiation between oxidative and fermentative carbohydrate metabolism.  
Perform an O.F. test (glucose) as described by MacFaddin (1980). Alternatively, inoculate tubes of GMD medium and interpret the results as described by Walters and Plumb (1978).
- d) Indole production  
Perform the test for indole production (MacFaddin 1980).
- e) Triple sugar iron agar  
The use of TSI agar and the interpretation of the results of this medium are described by MacFaddin (1980).
- f) Confirmatory slide agglutination test for A. salmonicida (Rabb et al. 1964) and Y. ruckeri.  
The agglutination test is performed by emulsifying a small amount of bacterial growth in saline (0.9% NaCl) on a clean glass slide. A loopful of antiserum is placed adjacent to the bacterial suspension and the two are mixed by gentle rocking and tilting of the slide. Prepare the appropriate positive and negative controls. A rapid macroscopic clumping of bacterial cells in the test mixture and the positive control (but not in the negative control) constitutes a positive agglutination test. Agglutination in the negative control invalidates the test. Many strains of A. salmonicida autoagglutinate. To prevent autoagglutination place the suspension in boiling water for 15 min prior to performing the slide agglutination test.

- g) Confirmatory tests for *R. salmoninarum*.  
If clinical signs or examinations of Gram-stained kidney smears suggest the presence of *R. salmoninarum*, retain a portion of the suspect kidney tissue and perform one of the following confirmatory tests. If such kidney tissue is unavailable, the presumptive evidence for *R. salmoninarum* noted earlier (IX B 2) constitutes a confirmed diagnosis.
- (1) Immunodiffusion test (Chen et al. 1974)  
Prepare immunodiffusion plates by adding ten mL of medium consisting of Noble agar (1.0%), NaCl (0.9%), and thimerosal (0.01%) to a 60-mm petri plate. Punch a pattern of six mm diameter wells, six of them peripheral to a center well, so that all wells are six mm apart. Place 0.1 mL specific antiserum in the center well. Load the peripheral wells separately with 0.1 mL of saline (negative control), a heavy *R. salmoninarum* cell suspension in saline (positive control), and a 50% kidney homogenate in saline (test sample). Arrange the loading so test samples are adjacent to positive controls. Incubate the plates in a moist chamber at 15°C for 48 h. A precipitin line of identity or partial identity between a test sample and a positive control constitutes a positive test.
  - (2) Direct fluorescent antibody test (DFAT) (Bullock et al. 1980)  
Prepare a kidney smear on a clean glass slide, allow to air dry and fix for 5-8 min in acetone at 20°C. Add 1-2 drops of the recommended optimal dilution of conjugated anti *R. salmoninarum* serum containing a 1:150-1:200 dilution of rhodamine counterstain (Difco) to the slide and allow to react for 5-8 min at 20-25°C. Rinse the slide and wash for two min in phosphate buffered saline (pH 7.2) and air dry. Add a drop of mounting fluid (pH 9.0) to the test area, add a coverslip and examine a minimum of 25 fields under oil immersion using a microscope equipped with an ultraviolet light source. A positive control should be prepared and stained in a similar manner. The presence of small fluorescing diplobacilli of typical size and shape constitutes a positive test.
  - (3) Indirect fluorescent antibody test (IFAT) (Bullock and Stuckey 1975)  
The indirect fluorescent antibody test can be used in place of the immunodiffusion test and DFAT to confirm the presence of *R. salmoninarum*. Consult the reference for the proper procedure.
- h) Biochemical confirmation of *A. salmonicida* and *Y. ruckeri*  
Confirmatory testing of presumptively identified isolates of *A. salmonicida* and *Y. ruckeri* is performed using conventional media (Difco) as described by Edwards and Ewing (1972) or the API-20E miniaturized diagnostic system (Analytab Products Inc., Plainview, NY, U.S.A.) and comparing the results with those obtained for known (positive control) cultures of the disease agent.

- (1) Streak bacterial growth on TSA and incubate at 20°C for 24-48 h to obtain a pure culture.
  - (2) Prepare and inoculate separately the suspect bacterial culture and the known culture as recommended (Difco, Analytab Products Inc.). When confirming presumptively identified isolates of Y. ruckeri, a saline suspension having a final turbidity equivalent to that of a #1 McFarland turbidity standard is recommended as the inoculum.
  - (3) Incubate the inoculated biochemical test media at 20°C for 24-72 hours. Add the necessary reagents and read the test results as recommended (Difco, Analytab Products Inc.).
  - (4) Compare the results obtained for the suspect isolate to those of the known bacterial culture(s). For interpretation of the biochemical profiles for A. salmonicida consult the papers by Paterson (1974) and Paterson et al. (1980); and for Y. ruckeri consult the paper by Stevenson and Daly (1982).
- i) Sensitivity to the vibriostatic agent 0/129 and the antibiotic novobiocin  
 The test is performed on TSA plates by applying a 0/129 disk and a novobiocin (5 µg) disk (Difco) to the medium that has been uniformly surface-seeded with the organism under test. After incubation at 20-22°C for 16-24 h, sensitive organisms show a clear zone around the disk. To prepare 0/129 disks, saturate Whatman antibiotic assay filter paper disks (6 mm) with 0.1% 0/129 in acetone. Drain off excess solution and dry the disks at 37°C. The vibriostatic agent 0/129 (2, 4-diamino-6, 7-diisopropylpteridine) is obtainable from British Drug Houses in Canada and from Calbiochem in the U.S.A. A control disk impregnated only with acetone should be included to preclude the possible inhibitory reaction due to acetone.
- j) Agar block motility test for myxobacteria
- (1) Excise a five mm square block of agar supporting a suspected myxobacterial colony. Place the block, colony side up, on a glass slide and cover gently with a cover slip.
  - (2) Examine the margin of the colony under high power for evidence of gliding or creeping motility.
- k) Confirmatory serological tests  
 Confirmation of identity of bacteria by serological tests such as slide, tube, or micro-well agglutination reactions should be carried out using standard procedures. Sera used should be standardized, preferably adsorbed and appropriate for the purpose. Positive and negative control organisms must be included.

## X PROCEDURES FOR THE DETECTION OF VIRUSES

### A. SCOPE

1. Any filterable agent in the fish samples that replicates intracellularly in any of the specified cell lines is certifiable whether or not it can be identified with presently available antisera and whether pathogenicity for salmonids exists or is unknown. Methodology is dependent upon detection of cytopathic effects (CPE) in susceptible cell cultures.
2. Any abnormal proliferative lesions (tumors) encountered should be processed by histological methods, and the results of the histopathologic evaluation reported.

### B. TISSUES TO BE ASSAYED

1. Sac fry: assay whole. When present, yolk sacs should first be removed and discarded.)
2. Fish averaging 2-4 cm in fork length: remove and discard heads, but retain gills; cut off the tails just posterior to the vents. Mince the remainder of the carcasses and assay.
3. Fish averaging 4-10 cm in fork length: excise the gills, then eviscerate and assay the combined viscera and gills. After removal of the gills, evisceration is readily accomplished by first cutting off the head, then slitting open the body cavity from the cut end to the vent, and finally cutting and scraping to remove the viscera (including kidney).
4. a) Fish averaging 10 cm or more in fork length: assay mixtures of kidney, spleen, pyloric caeca-pancreas, and gill. The appropriate relative volumes of these tissues should be 3:1:1:1, respectively. One part each of anterior, central, and posterior kidney must be represented in the sample.  
b) When the fish in this size category are brood fish, and reproductive fluids are to be used, as many as possible of the reproductive fluid samples must be from female fish. For maximum sensitivity, assay reproductive fluid samples individually.

### C. POOLING

Tissues from a maximum of five fish may be pooled to form one sample. But when preparing pooled samples, apparently healthy fish (or their tissues) must not be pooled with dead and moribund fish (or their tissues).

#### D. PREPARING INOCULA

Samples must be processed within 48 h of collection (see VIII A 1). Prior to and during processing, samples must be kept refrigerated or on ice, but not frozen.

1. Solid tissues: weigh and then homogenize the tissues in a minimal volume of balanced salt solution (BSS) such as Earle's or Hanks' BSS at pH 7.6-7.8. Two methods are available:

- a) Use a sterile Ten Broeck, or homogenizer designed to allow cooling on ice during homogenization to process small fish or small amounts of tissue. Care should be taken to prevent possible aerosol dissemination of virus.

- b) Use a sterile pre-chilled mortar and pestle to grind the tissues with a small quantity of sterile sand (80-120 mesh silica) until a smooth paste is formed.

Equipment used for homogenization of pooled samples within any given lot need not be sterilized between each use.

After tissues have been triturated, dilute each sample extract to a final concentration of 2% tissue suspension in BSS. Centrifuge the extracts at 2500 x g for 15 min at 4°C and aseptically filter the clarified supernatant through a 0.45 µm pore diameter membrane filter. To avoid any appreciable loss of virus by adsorption on the filter, collect the maximum possible volume of filtrate.

2. Fluid samples: dilute 1:2 with cold Eagle's minimum essential medium (MEM) at pH 7.2-7.6. Centrifuge at 2500 x g for 15 min at 4°C. Decontaminate as in X D 1.

#### E. ASSAY

1. Cell cultures

For virus assay use two of the following three continuous cell lines: rainbow trout gonad (RTG-2), chinook salmon embryo (CHSE-214), and fathead minnow (FHM). The FHM line must be used in IHN virus enzootic areas.

Biannually or before each inspection season all stock cell cultures should be tested and found to be free of mycoplasma. Each cell type must also be tested for susceptibility to viruses enzootic to the region (i.e. state, province or watershed). The donor cell cultures used to prepare monolayers for virus detection must be no older than two weeks.

Additional information on fish cell culture and virology can be found in McAllister (1979), Pilcher and Fryer (1980), Wolf and Quimby (1969) and Wolf (1970).

2. Either of the following procedures may now be used:

- a) Inoculation of preformed monolayers.

- (1) Prepare duplicate monolayers of each of two cell lines for each sample to be tested. Plastic multidishes

- suitable for tissue culture (1.5 to 2.0 cm diameter wells) may be used sealed or unsealed (with the appropriate organic buffer incorporated in the medium).
- (2) Use 1.0 mL per well of Eagle's MEM at pH 7.6-7.8 containing Earle's salts, glutamine, and 10% fetal bovine serum (FBS). Alternate antibacterial mixtures for use in the medium are 100 IU penicillin/mL and 100 µg streptomycin/mL, or 50 µg gentamycin/mL. The use of a fungistat (e.g., 25 IU nystatin/mL) is also permitted.
  - (3) Incubate the cell cultures at 15-20°C; the temperature depends upon when they will be required for the assay. At the time of inoculation the cell monolayers must be 70-90% confluent and not more than 48 h old.
  - (4) The growth medium must be removed and the monolayers washed with BSS, which is removed prior to inoculating 0.1 mL of the filter sterilized sample into each well.
  - (5) Incubate the inoculated cell cultures at 15°C for 60-90 min. Every 20 min gently rock cultures to uniformly spread the inocula. Add one mL of Eagle's MEM containing 2% FBS to each monolayer. Note: the new medium should be of the same composition as in E 2 a (2) except for decreased FBS and the final pH must be 7.6 to 7.8.
  - (6) Incubate cultures at 15°C.
- b) Simultaneously applied cells and test sample.
- (1) Place in tissue culture wells 1.0 mL of medium (see X E 2 a) containing sufficient cells to produce a 70-90% confluent monolayer on attachment. Duplicate cultures of two cell lines (see X E 1) are required per sample.
  - (2) Immediately add 0.1 mL of the filter sterilized sample to each culture.
  - (3) Incubate cultures at 15°C.

### 3. Controls

For each batch of donor cell cultures, duplicate negative controls must be run. Negative controls must consist of cultures inoculated according to the procedure used for the assays, except that sterile BSS must be used in place of the sample.

### 4. Procedures to be followed during incubation

- a) Inspect the cultures shortly after inoculation and after 24 h. Thereafter examine cultures at least every day to determine whether CPE has been produced.
- b) A sample is considered negative if there is no CPE in the cultures 14-21 days postinoculation.
- c) If CPE occurs in one or more of the cultures inoculated with the samples, the presence of a filterable, replicating agent must be verified. Filter (0.45 µm pore diameter) the culture fluid from the test well(s) showing CPE, dilute the filtrate

with BSS to  $10^{-1}$  and  $10^{-3}$ , and inoculate 0.1 mL of each of the dilutions into fresh duplicate cultures of the same cell line. If CPE is again observed, proceed with the serum neutralization test.

#### F. SERUM NEUTRALIZATION TEST

Presumptive identification of the agent producing CPE may be made on the basis of clinical evidence at sampling and on the type of CPE produced. Identification is accomplished by neutralization of the agent with specific antiserum. Failure to obtain any degree of neutralization using antisera prepared against known viruses will usually indicate the presence of a previously unrecognized virus or an atypical serotype of a known salmonid virus.

##### 1. Procedure

- a) Use a dilution of antiserum sufficient to neutralize an equal volume of a suspension containing  $10^2$ - $10^3$  TCID<sub>50</sub> per mL of the homologous virus.
- b) Filter the fluid from a culture showing CPE through a 0.45  $\mu$ m pore diameter membrane filter. Dilute the filtrate  $10^{-2}$  and  $10^{-6}$  with sterile BSS.
- c)
  - (1) Mix 0.3 mL of the diluted antiserum with 0.3 mL of each of the dilutions of sample.
  - (2) Mix 0.3 mL of normal serum with 0.3 mL of each of the dilutions of test sample.
  - (3) In the same manner, perform the serum neutralization test on the positive controls using homologous antiserum and normal serum.
- d) Incubate the reaction mixtures at 15°C for 30-60 min and then inoculate 0.2 mL of each mixture into duplicate cultures of the cell line in which the virus was isolated.
- e) Incubate the cultures at 15°C and observe for the production and inhibition of CPE. Inhibition of CPE by a particular antiserum, but not by normal serum, identifies the virus.

#### G. OTHER CONFIRMATORY METHODS

The approved methodology for detection of viruses is based upon isolation, followed by serological identification. Methods of confirming the identification of the cell culture isolates are not limited to the suggested serum neutralization test (X F). Other acceptable immunoserological tests may be used, including fluorescent-antibody microscopy, immunoperoxidase techniques, enzyme-linked immunosorbent assays, microtitration and microneutralization, plaque neutralization, complement fixation and immunoelectron microscopy.

## XI PROCEDURES FOR THE DETECTION OF CERTAIN PARASITES

### A. SCOPE

The absence of two myxosporean disease agents must be verified. These agents are Myxosoma cerebralis and Ceratomyxa shasta. Any other parasites detected which the Fish Health Official considers important should be reported.

### B. PROCEDURES FOR MYXOSOMA CEREBRALIS

1. Fish for this procedure must be at least 120 days old, and preferably fresh. Frozen, but not formalin-preserved specimens may also be used. All glassware and equipment used in processing samples must be carefully cleaned to avoid carry-over of spores.
2. Decapitate fish and deflesh heads after heating in water at 45-50°C until the brain has coagulated. Remove the brain (and any attached spinal cord) intact. Discard to avoid possible contamination of skull material with spores of the related parasite Myxobolus neurobius that might be mistaken for M. cerebralis.
3. Either of the following procedures may be now used:
  - a) Digestion method<sup>4</sup>
    - (1) Pool up to approximately 100 g of the resulting skull material and macerate finely. Small initial fragment size will facilitate complete and rapid digestion.
    - (2) Place in a beaker and add 25 mL of freshly prepared pepsin digest solution (1.0 g of pepsin dissolved in 100 mL of 0.5% HCl) to each gram of macerated material.
    - (3) Mix well, let settle for 2 min and examine a sample from the surface of the supernatant at 400-450X magnification for typical spores using phase contrast microscopy.
    - (4) If spores are not detected, incubate the mixture at 35-40°C for 1-1.5 h with gentle agitation. A cloudy, greyish suspension, free of large particles, should be present at the end of the digest period.
    - (5) Place 50 mL of the suspension into conical, screw cap centrifuge tubes designed to hold 50 mL. Centrifuge at 1200 x g for 15 min at room temperature. Discard the supernatant and resuspend the pellet in 1.0 mL of

---

<sup>4</sup>The procedure outlined above is a modification of the pepsin-trypsin digest method of Markiw and Wolf (1974a and 1974b).

distilled water. The contents of up to five tubes can be combined into one sample. Examine microscopically for spores as before.

- (6) If no spores are detected, bring the volume of each sample to approximately six mL with distilled water. Layer the contents of each tube onto three mL of 55% aqueous glucose solution contained in a 12 mL conical centrifuge tube. Centrifuge at 1200 x g for 30 min at room temperature.
  - (7) Withdraw pelleted material with a Pasteur pipette and examine at least 25 microscope fields per sample as given in XI B 3 a (3). Observation of spores at any stage in the procedure constitutes a positive result.
- b) Plankton centrifuge method (O'Grodnick 1975 and Prasher et al. 1971)
- (1) Pool up to 100 g of skull material and macerate in a blender for five min with distilled water. Up to 200 mL of water per 100 g of skull material can be used.
  - (2) Remove the macerated material and vacuum filter through a fine (0.5-1.0 mm) wire mesh. If the filter clogs, rinse with distilled water to clear it allowing the rinse water to mix with the filtrate. Coarse material such as large bone chips are removed by this procedure and can be discarded.
  - (3) Place the filtrate in a separatory funnel located to discharge into a plankton centrifuge. Run the centrifuge at high speed while adding the filtrate in slowly.
  - (4) Centrifuge until all water has been removed. Scrape the residue from the wall of the centrifuge and place in a small bottle. Add five volumes of distilled water. Do not dilute beyond a total of 30 mL. Cap and shake the bottle until the material is uniformly suspended.
  - (5) Place a drop of the suspension on a slide or hemocytometer if quantitation is required.
  - (6) Examine at least 25 fields of the slide at 450X magnification for the presence of M. cerebralis spores using phase contrast microscopy.

4. Confirmatory identification of M. cerebralis spores must be based on the morphological characters given by Lom and Hoffman (1971).

#### C. PROCEDURES FOR CERATOMYXA SHASTA

1. Fish for the following procedures must be at least 120 days old; either fresh (preferably) or frozen specimens may be used.
2. The preferred organs for examination for C. shasta spores are the intestine and gall bladder. Peritoneal fluid may also contain

spores. Microscope slides of tissue, fluid and purulent material can be prepared and examined as required by either of two methods:

a) Wet mounts

Prepare wet mounts by gently mixing sufficient material in one or two drops of saline (0.9% NaCl) on a standard microscope slide to give a reasonably dilute suspension, cover with a cover slip and examine at 400-450X magnification, using phase contrast microscopy.

b) Dried smears

Smear material to be examined on a standard microscope slide, allow to air dry, stain 30-60 sec with Loeffler's methylene blue (dissolve 0.3 g of methylene blue chloride in 30 mL of 95% ethanol and add 100 mL of 0.1% aqueous KOH), rinse with water and air dry. Add a drop of immersion oil or water to each smear, cover with a cover slip and examine at 400-450X magnification for the presence of spores. The polar capsules and extended polar filaments stain an intense blue and the sporoplasm stains a pale blue.

3. For fish less than 7.5 cm in fork length, material for examination can be expressed from the intestines. In addition, if nodular lesions are present in any tissue, or if ascitic fluid is evident, prepare smears or wet mounts of this tissue or fluid.
4. For fish 7.5 cm or more in fork length, open the body cavity and collect material lightly scraped from the interior wall of the upper intestine or gall bladder. If nodular lesions are observed on or in any tissue, especially the pyloric ceca, or if any abnormal accumulations of fluid are found, they must be examined for spores of C. shasta.
5. For each smear or wet mount at least 25 microscope fields must be examined for spores of C. shasta.
6. The identification of C. shasta spores must be based on the diagnostic characters given by Johnson et al. (1979). Pre-spore stages of C. shasta may be found without accompanying spores; by themselves they are not diagnostic for the organism and their presence indicates that further samples and search for the characteristic spores should be made (Noble 1950; Yamamoto and Sanders 1979).

Additional information on identification of fish disease agents can be found in McDaniel (1979).

## XII EGG DISINFECTION PROCEDURES

Salmonid eggs are safely disinfected as green eggs following fertilization and water hardening, or as early eyed eggs. The following is a suggested procedure for egg disinfection, utilizing iodophors. Iodophors for disinfection are usually povidone or polyalcoholic complexes of iodine in which the solubilized iodine confers its broad spectrum germicidal activity but is not as corrosive or irritating as in its elemental form. A number of topical disinfectants<sup>5</sup> of this type are available commercially in North America; among these are Ovadine®, Bridine®, Betadine®, Actomar K30®, Wescodyne® and Argentyne®. Most contain a 1-2% active iodine concentration.

### A. PREPARATION OF THE DISINFECTANT

1. Dilute the stock iodine-based disinfectant to give a solution containing 100 parts per million (ppm) of active iodine. The disinfectant must be prepared in water with a low organic content to minimize loss of the free iodine. Use a plastic, glass, stainless steel or fiberglass tank for preparing the holding the solution.
2. Check the pH of the diluted disinfectant and, if necessary, adjust to 6.5-7.5 using 8% aqueous sodium bicarbonate (baking soda).

### B. DISINFECTION PROCEDURE

1. Use a fresh solution of diluted disinfectant.
2. To avoid temperature shock, adjust the disinfectant solution to the same temperature as the subsequent egg incubation temperature.
3. In the case of freshly fertilized eggs, allow eggs to water harden one hour before disinfection.
4. Immerse water hardened green eggs or early eyed eggs in the disinfectant for ten min.
5. Treat approximately 2000 eggs per liter before discarding the disinfectant.
6. Rinse eggs thoroughly in uncontaminated water after disinfection.

---

<sup>5</sup>The products specified have proven satisfactory for the purposes indicated; this, however, does not imply that other products may not be equally satisfactory.

7. Arrange the egg handling program to ensure that disinfected eggs do not have subsequent contact with contaminated equipment, water or personnel.

Diluted iodophors can also be used to disinfect work surfaces, utensils, nets and other equipment used during the egg taking process, but rinse thoroughly in clean, uncontaminated water following the disinfection.

## REFERENCES

- AMANDI, A., S.F. HIU, J.S. ROHOVEC and J.L. FRYER. 1982. Isolation and characterization of Edwardsiella tarda from fall chinook salmon (Oncorhynchus tshawytscha). Appl. Envir. Microbiol. 43: 1380-1384.
- ANACKER, R.L., and E.J. ORDAL. 1959. Studies on the myxobacterium Chondrococcus columnaris. 1. Seriological typing. J. Bacteriol. 78: 25-32.
- BULLOCK, G.L., B.R. GRIFFIN and H.M. STUCKEY. 1980. Detection of Corynebacterium salmoninus by direct fluorescent antibody test. Can. J. Fish. Aquat. Sci. 37: 719-721.
- BULLOCK, G.L. and H.M. STUCKEY. 1975. Fluorescent antibody identification and detection of the Corynebacterium causing kidney disease of salmonids. J. Fish. Res. Board Can. 32: 2224-2227.
- CHEN, P.K., G.L. BULLOCK, H.M. STUCKEY, and A.C. BULLOCK. 1974. Serological diagnosis of corynebacterial kidney disease of salmonids. J. Fish. Res. Board Can. 31: 1939-1940.
- EDWARDS, P.R. and W.H. EWING. 1972. Identification of Enterobacteriaceae, 3rd ed. Burgess Publishing Co., Minneapolis, MN. 362 p.
- EVELYN, T.P.T. 1971. An aberrant strain of the bacterial fish pathogen Aeromonas salmonicida isolated from a marine host, the sablefish (Anoplopoma fimbria) and from two species of cultured Pacific salmon. J. Fish. Res. Board Can. 28: 1629-1634.
- FIJAN, N.N. 1969. Antibiotic additives for the isolation of Chondrococcus columnaris from fish. Appl. Micro. 17: 333-334.
- GRIFFIN, P.J., S.F. SNIESZKO, and S.B. FRIDDLE. 1952. A more comprehensive description of Bacterium salmonicida. Trans. Am. Fish. Soc. 82: 129-138.
- JOHNSON, K.A., J.E. SANDERS, and J.L. FRYER. 1979. Ceratomyxa shasta in salmonids. U.S. Fish and Wild. Serv., Fish Dis. Leaflet. No. 58, Washington, DC. 11 p.
- LOM, J. and G.L. HOFFMAN. 1971. Morphology of the spores of Myxosoma cerebralis and M. cartilaginis (Hoffman, Putz, and Dunbar 1965). J. Parasitol. 57(6): 1302-1308.
- MacFADDIN, J.F. 1980. Biochemical tests for the identification of medical bacteria, 2nd ed. Williams and Wilkins, Baltimore, MD. 527 p.
- MARKIW, M.E., and K. WOLF. 1974a. Myxosoma cerebralis: Isolation and concentration from fish skeletal elements - sequential enzymatic digestions and purification by differential centrifugation. J. Fish. Res. Board Can. 31: 15-20.
- 1974b. Myxosoma cerebralis: Comparative sensitivity of spore detection methods. J. Fish. Res. Board Can. 31: 1597-1600.
- McALLISTER, P.E. 1979. Fish viruses and viral infections, In: Comprehensive Virology, Vol. 14, Fraenkel-Conrat, H. and Wagner, R.R. eds., Plenum Press, New York and London. 401-470.
- McDANIEL, D. (ed.). 1979. Fish health bluebook: procedures for the detection and identification of certain fish pathogens. Fish. Health Sec., Am. Fish. Soc., Bethesda, MD. 118 p.

- NOBLE, E.R. 1950. On a myxosporidian (protozoan) parasite of California trout. *J. Parasitol* 36: 457-460.
- O'GRODNICK, J.J. 1975. Whirling disease Myxosoma cerebralis: Spore concentration using the continuous plankton centrifuge. *J. Wildl. Dis.* 11: 54-57.
- OSSIANDER, F.J., and G. WEDEMEYER. 1973. Computer program for sample sizes required to determine disease incidence in fish populations. *J. Fish. Res. Board Can.* 30: 1383-1384.
- PATERSON, W.D. 1974. Biochemical and serological differentiation of several pigment producing aeromonads. *J. Fish. Res. Board Can.* 31: 1259-1261.
- PATERSON, W.D., D. DOUEY, and D. DESAUTELS. 1980. Isolation and identification of an atypical Aeromonas salmonicida strain causing epizootic losses among Atlantic salmon (Salmo salar) reared in a Nova Scotian hatchery. *Can. J. Fish. Aquat. Sci.* 37: 2236-2241.
- PILCHER, K.S., and J.L. FRYER. 1980. The viral diseases of fish: A review through 1978 Part 1: Diseases of proven viral etiology. *CRC Critical Reviews in Microb.* 7(4): 287-363.
- PRASHER, J.B., W.M. TIDD, and R.A. TUBB. 1971. Techniques for extracting and quantitatively studying the spore stage of the protozoan parasite Myxosoma cerebralis. *Prog. Fish-Cult.* 33: 193-196.
- RABB, L., J.W. CORNICK, and L.A. McDERMOTT. 1964. A microscopic slide agglutination test for the presumptive diagnosis of furunculosis in fish. *Prog. Fish-Cult.* 26: 118-119.
- SANDERS, J.E. and J.L. FRYER. 1980. Renibacterium salmoninarum gen. nov., sp. nov., the causative agent of bacterial kidney disease. *Int. J. Syst. Bacteriol.* 30(2): 496-502.
- SHIEH, H.S. 1980. Studies on the nutrition of a fish pathogen, Flexibacter columnaris. *Microbios Letters*, 13: 129-133.
- STEVENSON, R.M.W. and J.G. DALY. 1982. Biochemical and serological characteristics of Ontario isolates of Yersinia ruckeri. *Can. J. Fish. Aquat. Sci.* 39: 870-876.
- TAYLOR, E.L., S.J. COLI, and D.R. JUNELL. 1973. Attempts to control whirling disease by continuous drug feeding. *J. Wildl. Dis.* 9: 320-325.
- WOLF, K. 1970. Guidelines for virological examination of fishes, In: A symposium on diseases of fish and shellfish. *Am. Fish. Soc., Spec. Publ. No. 5*, Washington, DC. 327-340.
- WOLF, K. and M.C. QUIMBY. 1969. Fish cell and tissue culture, In: *Fish Physiol.* Vol. 3, Hoar, W.S. and Randall, D.J., eds. Academic Press, New York, NY. 253-305.
- YAMAMOTO, T. and J.E. SANDERS. 1979. Light and electron microscopic observations of sporogenesis in the myxosporidium, Ceratomyxa shasta (Noble 1950). *J. Fish Dis.* 2: 411-428.

## APPENDIX 1 NATIONAL REGISTRY OF FISH DISEASES

The National Registry of Fish Diseases is a data center for the documentation and dissemination of information pertaining to fish diseases in Canada. Fish Health Laboratory Reports of all disease examinations in Canada, not only those specifically for certification purposes but also those from routine periodic sampling programs should be submitted routinely to the Registry by Fish Health Officials.

The National Registry of Fish Diseases will:

1. maintain close watch on geographic distribution and prevalence of fish diseases in Canada and identify emerging disease problems;
2. serve as a coordinating center in the event of national fish health emergencies;
3. provide periodic reports, analyses, and assessments on the state of fish health in Canada;
4. provide health histories of sources of live fish and eggs;
5. maintain and provide lists of certified Canadian and foreign fish production facilities;
6. maintain and provide lists of Local Fish Health Officers and Fish Health Officials approved by the Minister; and
7. maintain and provide lists of Canadian government approved, certifying officials in foreign countries involved in export of fish and eggs to Canada.

The address for the Registry is:

The National Registry of Fish Diseases,  
Department of Fisheries and Oceans,  
Ottawa, Ontario, Canada  
K1A 0E6 (Tel. 613-990-0276; Telex 0534228)

## APPENDIX 2 REGIONAL ADMINISTRATIVE AUTHORITIES

The Fish Health Protection Regulations are administered and applied by various provincial and regional authorities. In some areas provincial authorities are involved while other areas are a federal responsibility. The addresses follow:

### 1. Provincial fisheries authorities

#### Alberta

Fish and Wildlife Division,  
Department of Energy and Natural Resources,  
9915 - 108 Street,  
Edmonton, Alta.  
T5K 2C9 (Tel. 403-427-6730)

#### British Columbia

Fish and Wildlife Branch,  
Ministry of Environment,  
780 Blanshard Street,  
Victoria, B.C.  
V8V 1X5 (Tel. 604-387-1961)

#### Manitoba

Fisheries Enhancement Services,  
Fisheries Branch,  
Department of Natural Resources,  
1495 St. James Street,  
Winnipeg, Man.  
R3G 0W9 (Tel. 204-944-7792; Telex 07-587740)

#### Ontario

Fisheries Branch,  
Ontario Ministry of Natural Resources,  
Queen's Park,  
99 Wellesley St. West,  
Toronto, Ont.  
M7A 1W3 (Tel. 416-965-7886; Telex 06-219701)

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#### Saskatchewan

Superintendent of Fisheries,  
Fisheries Branch,  
Department of Parks and Renewable Resources,  
Box 3003,  
Prince Albert, Sask.  
S6V 6G1 (Tel. 306-922-9887)

2. Regional federal fisheries authorities

Manitoba

Inspection and Technology Branch,  
Department of Fisheries and Oceans,  
153 Lombard Ave.,  
Winnipeg, Man.  
R3B 0T3 (Tel. 204-949-4060; Telex 075-7419)

Newfoundland

Newfoundland Regional Office,  
Fisheries Research Branch,  
Department of Fisheries and Oceans,  
P.O. Box 5667,  
St. John's, Nfld.  
A1C 5X1 (Tel. 709-772-2194; Telex 016-4698)

Northwest Territories

Department of Fisheries and Oceans,  
Box 1008, Hay River,  
N.W.T. (Tel. 403-874-2331; Telex 034-4324)

New Brunswick, Nova Scotia, Prince Edward Island

Scotia-Fundy Regional Office,  
Department of Fisheries and Oceans,  
P.O. Box 550,  
Halifax, N.S.  
B3J 2S7 (Tel. 902-426-2581; Telex 019-21819)

Gulf Regional Office,  
Department of Fisheries and Oceans,  
P.O. Box 5030,  
Moncton, N.B.  
E1C 9B6 (Tel. 506-758-9044; Telex 0142-607)

Ontario

Fishing and Industry Services,  
Department of Fisheries and Oceans,  
590 Keele Street, Room 410,  
Toronto, Ont.  
M6N 3E3 (Tel. 416-763-1161; Telex 06-22790)

Québec

Inspection and Technology Branch,  
Department of Fisheries and Oceans,  
1001 Pierre Dupuy,  
Captain Bernier Laboratory,  
Longueuil, Que.  
J4K 1A1 (Tel. 514-283-4967)

Saskatchewan

District Manager,  
Department of Fisheries and Oceans,  
Immigration Building,  
Prince Albert, Sask.  
S6V 1E8 (Tel. 306-764-6424; Telex 074-29217)

Yukon and marine areas of British Columbia

Pacific Biological Station,  
Department of Fisheries and Oceans,  
Nanaimo, B.C.  
V9R 5K6 (Tel. 604-756-7000; Telex 044-6128)

### APPENDIX 3 QUALIFICATIONS OF FISH HEALTH OFFICIALS

A Fish Health Official carries the burden of a great deal of responsibility and therefore must be knowledgeable in the field of fish health and thoroughly familiar with the methodology for diagnosing fish diseases and their causative infectious agents. The following criteria represent the minimum requirements for acceptance by the Canadian government as a Fish Health Official.

#### A. EDUCATION AND EXPERIENCE

1. Bachelor's degree or equivalent in any field of fisheries or biological science and two years' experience in the diagnosis of fish diseases, including experience in bacteriological, virological, and parasitological techniques.
2. Veterinary, master's, or doctoral degree or equivalent with evidence of training in bacteriology, virology, and parasitology and one year's experience in the diagnosis of fish diseases.

#### B. LABORATORY FACILITIES

Laboratory facilities adequate for microbiological and parasitological evaluation of samples must be available to the proposed Fish Health Official.

#### C. APPLICATION PROCEDURE

To be considered for approval as a Fish Health Official, a candidate must submit a curriculum vitae, three specimen signatures, and a description of available laboratory facilities and equipment to the National Registry of Fish Diseases.

APPENDIX 4

Fish Health Certificate

Government of Canada  
Fisheries and Oceans

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Certificate Expiry Date  
D M Y

#### CERTIFICATION OF SOURCE

Name of source	Address	Telephone No.
----------------	---------	---------------

I \_\_\_\_\_, (Print name), authorized as a Fish Health Official by the Minister of Fisheries and Oceans of Canada; certify that the source indicated above was inspected by methods approved by the Minister, and that no evidence was found for the diseases or disease agents for  livefish (Schedule II) or for  dead fish (Schedule III) as required by the Fish Health Protection Regulations made under the Fisheries Act of Canada.

The following notifiable disease agents (Schedule IV) and/or other fish pathogens were found: \_\_\_\_\_

Dates of last four previous satisfactory inspections	D M Y	D M Y	D M Y	D M Y
--	-------	-------	-------	-------

Date of issue	Signature and address of Fish Health Official	Telephone No.
---------------	---	---------------

#### CONSIGNMENT

Departing city and country	Carrier	Bill of lading No.	Date D M Y
----------------------------	---------	--------------------	---------------

Anticipated port of arrival in Canada (City and Province)	Date D M Y
---	---------------

Date	Name and address of person receiving shipment	Telephone No.
------	---	---------------

#### EXPORTER'S DECLARATION

I \_\_\_\_\_, (Print name),  owner  manager of the above noted facility which was last inspected on \_\_\_\_\_ Date, declare that to my knowledge no disease agents listed in Schedule II of the Fish Health Protection Regulations (FHPR) have been detected according to the procedures outlined in the FHPR Manual of Compliance in this facility since the last FHPR inspection, that no introduction of fish or fish eggs from an uncertified source has been made into the facility, that the shipment described above will be derived solely from this facility, and that eggs in the shipment will be surface disinfected prior to leaving the source.

I \_\_\_\_\_, (Print name), consignor of eggs taken from wild spawners declare that these eggs will be surface disinfected and that they derive solely from the above inspected source.

This shipment consists of:

_____ kg	<input type="checkbox"/> Live	<input type="checkbox"/> Eggs	_____ Species
_____ Number	<input type="checkbox"/> Dead	<input type="checkbox"/> Fish	_____ Species

Date	Signature and address of owner, manager or consignor	Telephone No.
------	--	---------------

**APPENDIX 5**

**Fish Health Laboratory Report**

Inspection Date	D	M	Y
FHPR Schedule			
Investigator			
Laboratory			

**FISH HEALTH LABORATORY REPORT**

A. Location Data	1. Name and address of facility		2. Type of facility			
			3. Natural water body (specify)	4. Province		
B. Fish Data	Species	Fish/Eggs	Age/Stage	Number held	Number examined	
Stocks introduced since last FHPR inspection						
	Species	Age/Stage	Number	Source		
C. Water Data	1. <input type="checkbox"/> Municipal water	2. <input type="checkbox"/> Lake	3. <input type="checkbox"/> Stream	4. <input type="checkbox"/> Enclosed spring		
	5. <input type="checkbox"/> Open spring	6. <input type="checkbox"/> Well	7. <input type="checkbox"/> Deepwell	8. <input type="checkbox"/> Brackish water		
	9. <input type="checkbox"/> Marine	10. <input type="checkbox"/> Influent disinfection	11. <input type="checkbox"/> Effluent disinfection			
	What does effluent water drain into?			Name of watershed		
D. Disease Data	1. <input type="checkbox"/> No pathogens detected		2. <input type="checkbox"/> Schedule II pathogens detected			
	Indicate prevalence by denoting: <input type="checkbox"/> carriers; <input type="checkbox"/> clinical signs or <input type="checkbox"/> epizootic in boxes below.					
	<input type="checkbox"/> <i>A. salmonicida</i> (furunculosis)	<input type="checkbox"/> VHSV	<input type="checkbox"/> <i>M. cerebralis</i>			
	<input type="checkbox"/> KD bacterium <i>R. salmoninarum</i>	<input type="checkbox"/> IHNV	<input type="checkbox"/> <i>C. shasta</i>			
<input type="checkbox"/> RM bacterium <i>Y. ruckeri</i>	<input type="checkbox"/> IPNV	<input type="checkbox"/> Other replicating agents				
3. Other pathogens or tumors detected (specify)						
_____						
Additional comments (e.g. facility disinfection and chemicals used; egg disinfection and chemicals used; drugs/vaccines used)						
_____						
_____						
Date of issue			Signature of Fish Health Official			

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#### Saskatchewan

Superintendent of Fisheries  
Fisheries Branch  
Department of Parks and Renewable Resources  
Box 3003  
Prince Albert (Sask.)  
S6V 6G1 (Tél. : 306-922-9887)

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- NOBLE, E.R. 1950. On a myxosporidian (protozoan) parasite of California trout. *J. Parasitol* 36 : 457-460.
- O'GRODNICK, J.J. 1975. Whirling disease Myxosoma cerebralis : Spore concentration using the continuous plankton centrifuge. *J. Wildl. Dis.* 11 : 54-57.
- OSSIANDER, F.J., and G. WEDEMEYER. 1973. Computer program for sample sizes required to determine disease incidence in fish populations. *J. Fish. Res. Board Can.* 30 : 1383-1384.
- PATERSON, W.D. 1974. Biochemical and serological differentiation of several pigment producing aeromonads. *J. Fish. Res. Board Can.* 31 : 1259-1261.
- PATERSON, W.D., D. DOUEY, and D. DESAUTELS. 1980. Isolation and identification of an atypical Aeromonas salmonicida strain causing epizootic losses among Atlantic salmon (Salmo salar) reared in a Nova Scotian hatchery. *Can. J. Fish. Aquat. Sci.* 37 : 2236-2241.
- PILCHER, K.S., and J.L. FRYER. 1980. The viral diseases of fish : A review through 1978 Part 1 : Diseases of proven viral etiology. *CRC Critical Reviews in Microb.* 7(4) : 287-363.
- PRASHER, J.B., W.M. TIDD, and R.A. TUBB. 1971. Techniques for extracting and quantitatively studying the spore stage of the protozoan parasite Myxosoma cerebralis. *Prog. Fish-Cult.* 33 : 193-196.
- RABB, L., J.W. CORNICK, and L.A. McDERMOTT. 1964. A microscopic slide agglutination test for the presumptive diagnosis of furunculosis in fish. *Prog. Fish-Cult.* 26 : 118-119.
- SANDERS, J.E. and J.L. FRYER. 1980. Renibacterium salmoninarum gen. nov., sp. nov., the causative agent of bacterial kidney disease. *Int. J. Syst. Bacteriol.* 30(2) : 496-502.
- SHIEH, H.S. 1980. Studies on the nutrition of a fish pathogen, Flexibacter columnaris. *Microbios Letters*, 13 : 129-133.
- STEVENSON, R.M.W. and J.G. DALY. 1982. Biochemical and serological characteristics of Ontario isolates of Yersinia ruckeri. *Can. J. Fish. Aquat. Sci.* 39 : 870-876.
- TAYLOR, E.L., S.J. COLI, and D.R. JUNELL. 1973. Attempts to control whirling disease by continuous drug feeding. *J. Wildl. Dis.* 9 : 320-325.
- WOLF, K. 1970. Guidelines for virological examination of fishes, In : A symposium on diseases of fish and shellfish. *Am. Fish. Soc., Spec. Publ. No. 5*, Washington, DC. 327-340.
- WOLF, K. and M.C. QUIMBY. 1969. Fish cell and tissue culture, In : *Fish Physiol. Vol. 3*, Hoar, W.S. and Randall, D.J., eds. Academic Press, New York, NY. 253-305.
- YAMAMOTO, T. and J.E. SANDERS. 1979. Light and electron microscopic observations of sporogenesis in the myxosporidium, Ceratomyxa shasta (Noble 1950). *J. Fish Dis.* 2 : 411-428.

## BIBLIOGRAPHIE

- AMANDI, A., S.F. HIU, J.S. ROHOVEC and J.L. FRYER. 1982. Isolation and characterization of Edwardsiella tarda from fall chinook salmon (Oncorhynchus tshawytscha). Appl. Envir. Microbiol. 43 : 1380-1384.
- ANACKER, R.L., and E.J. ORDAL. 1959. Studies on the myxobacterium Chondrococcus columnaris. 1. Seriological typing. J. Bacteriol. 78 : 25-32.
- BULLOCK, G.L., B.R. GRIFFIN and H.M. STUCKEY. 1980. Detection of Corynebacterium salmoninus by direct fluorescent antibody test. Can. J. Fish. Aquat. Sci. 37 : 719-721.
- BULLOCK, G.L. and H.M. STUCKEY. 1975. Fluorescent antibody identification and detection of the Corynebacterium causing kidney disease of salmonids. J. Fish. Res. Board Can. 32 : 2224-2227.
- CHEN, P.K., G.L. BULLOCK, H.M. STUCKEY, and A.C. BULLOCK. 1974. Serological diagnosis of corynebacterial kidney disease of salmonids. J. Fish. Res. Board Can. 31 : 1939-1940.
- EDWARDS, P.R. and W.H. EWING. 1972. Identification of Enterobacteriaceae, 3rd ed. Burgess Publishing Co., Minneapolis, MN. 362 p.
- EVELYN, T.P.T. 1971. An aberrant strain of the bacterial fish pathogen Aeromonas salmonicida isolated from a marine host, the sablefish (Anoplopoma fimbria) and from two species of cultured Pacific salmon. J. Fish. Res. Board Can. 28 : 1629-1634.
- FIJAN, N.N. 1969. Antibiotic additives for the isolation of Chondrococcus columnaris from fish. Appl. Micro. 17 : 333-334.
- GRIFFIN, P.J., S.F. SNIESZKO, and S.B. FRIDDLE. 1952. A more comprehensive description of Bacterium salmonicida. Trans. Am. Fish. Soc. 82 : 129-138.
- JOHNSON, K.A., J.E. SANDERS, and J.L. FRYER. 1979. Ceratomyxa shasta in salmonids. U.S. Fish and Wild. Serv., Fish Dis. Leaflet No. 58, Washington, DC. 11 p.
- LOM, J. and G.L. HOFFMAN. 1971. Morphology of the spores of Myxosoma cerebralis and M. cartilaginis (Hoffman, Putz, and Dunbar 1965). J. Parasitol. 57(6) : 1302-1308.
- MacFADDIN, J.F. 1980. Biochemical tests for the identification of medical bacteria, 2nd ed. Williams and Wilkins, Baltimore, MD. 527 p.
- MARKIW, M.E., and K. WOLF. 1974a. Myxosoma cerebralis : Isolation and concentration from fish skeletal elements - sequential enzymatic digestions and purification by differential centrifugation. J. Fish. Res. Board Can. 31 : 15-20.
- 1974b. Myxosoma cerebralis : Comparative sensitivity of spore detection methods. J. Fish. Res. Board Can. 31 : 1597-1600.
- McALLISTER, P.E. 1979. Fish viruses and viral infections, In : Comprehensive Virology, Vol. 14, Fraenkel-Conrat, H. and Wagner, R.R. eds., Plenum Press, New York and London. 401-470.
- MCDANIEL, D. (ed.). 1979. Fish health bluebook : procedures for the detection and identification of certain fish pathogens. Fish. Health Sec., Am. Fish. Soc., Bethesda, MD. 118 p.

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OL  
 626  
 C 314  
 NO. 31R  
 (1997)

**Regulations Amending the Fish Health Protection  
 Regulations**

*Statutory Authority*

*Fisheries Act*

*Sponsoring Department*

Department of Fisheries and Oceans

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*Ministère responsable*

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**REGULATORY IMPACT  
 ANALYSIS STATEMENT**

*Description*

The purpose of this amendment is to provide local fish health officers with more flexibility to issue import permits under the *Fish Health Protection Regulations*.

Under the current Regulations, live cultured eggs and fish and eggs of wild fish cannot be imported into Canada or transferred between provinces if any disease or disease agent set out in Schedule II to the Regulations is detected in the source facility.

If a receiving facility has the same diseases or disease agents as the source facility, there is no risk that a shipment of eggs or fish will result in the introduction of a new disease or disease agent to a province.

Hence the Regulations are being amended to allow local fish health officers to approve the importation of eggs or fish from a facility where a disease or disease agent of concern is present in both the source and receiving facilities. The guiding principle for local fish health officers will be that the importation or inter-provincial shipment of live eggs or fish will be allowed if the shipment will not result in the introduction of a disease or disease agent not already known to occur in a province. The same principles will apply to importation of dead un-eviscerated fish.

In the absence of data, a local fish health officer may designate facilities or areas in a province free of selected pathogens, even though the pathogens have been detected in other parts of a province. Hence, local fish health officers may reject applications to import eggs or fish, from a facility that tested positive for a selected pathogen, into a facility or area designated free of the specific pathogen.

Under the current Regulations, sources of salmonid eggs must be inspected for, and demonstrated free of, all diseases or disease agents set out in Schedule II. As a result of these amendments, local fish health officers would be able to approve the importation or inter-provincial transfer of disinfected salmonid

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eggs<sup>1</sup> from sources that have only been inspected for viruses (i.e. the sources have no history of testing for bacteria and parasites set out in Schedule II). The reason for this change is that there is minimal risk that bacterial and parasite disease agents listed in Schedule II will be carried on the outer membrane (surface) of disinfected eggs.

#### Alternatives

One alternative is to delay implementing this amendment until other major changes to the *Fish Health Protection Regulations* have been completed. These changes are expected to be ready for implementation in 1998. This would be unnecessarily restrictive on the aquaculture industry, because the amendment can be implemented immediately without increasing the risk of introducing or spreading disease agents.

#### Benefits and Costs

##### Benefits

Allowing importation or interprovincial transfers of live eggs and fish between source and recipient facilities that have the same fish disease profiles will increase the flexibility of aquaculturists to choose suitable sources of seedstock, and could increase the numbers of eggs and fish imported annually or shipped interprovincially. While these increases could result in increased requests to local fish health officers for import permits, administration costs should not increase substantially.

The impact of only requiring testing of sources of eggs for viruses will be two-fold. Firstly, facilities where bacterial or parasite disease agents of concern were previously detected (and hence they could not be "certified" under the current Regulations) could now become sources of eggs for export to Canada or for shipment between provinces. This could create new business opportunities for some private sector salmonid aquaculturists. Secondly, the cost of inspecting sources of eggs prior to the issue of fish health certificates will be substantially reduced (approximately 50 percent), because there will no longer be a requirement to test for bacteria and parasites. This is especially significant in Canada, where the private sector will soon be required to pay the costs for inspections leading to the issue of Fish Health Certificates. Canada does not pay for inspections of sources that export products into Canada.

##### Costs

There could be an increase in the number of facilities requesting inspection and disease testing for health certification. This would result in increased costs for health certification. However, the Department of Fisheries and Oceans (DFO) which now pays all the costs for health certification under the Regulations, intends to begin recovering these costs. So costs to the federal government should be reduced, not increased.

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<sup>1</sup> There are no chemicals registered in Canada for disinfection of salmonid eggs. However, the Department of Fisheries and Oceans is working with Health Canada and Agriculture and Agri-Food Canada to address this need.

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Costs for health certification of commercial egg suppliers should be reduced if eggs are disinfected, because facilities only have to be tested for viruses (not bacteria or parasites).

#### Consultation

There has been extensive consultation with government agencies, industry, and interested associations on proposals to change the Regulations. The two principles proposed (allowing transfers between facilities having equivalent health status; and reduced inspection requirements for diseases and disease agents in sources of eggs) are widely accepted by governments and industry.

Industry requested that an appeal procedure be set up to ensure that there is consistent implementation of the amendments nationally. A standard appeal procedure will be developed in consultation with fish health administrators and industry, and incorporated into the Manual of Compliance to the Regulations.

Concerns were expressed in some regions that the amendment could result in the transfer into a province of new strains of viruses or antibiotic resistant strains of bacteria. These concerns would be better addressed through development of regional or provincial policies designed to address specific local issues, rather than through national legislation.

Some government agencies and industry representatives recommended that there be a standard procedure, based on risk assessment principles, for identifying specific pathogen-free zones in a province. A standard procedure will be developed in consultation with fish health administrators and industry.

DFO will notify local fish health officers, fish health officials, aquaculture associations and individual producers by correspondence, giving details of the amendment and compliance procedures. Implementation of the amendment will be reviewed on a regular basis by local fish health officers from across the country. Information on Compliance Procedures may be obtained from the National Registry of Fish Diseases, Department of Fisheries and Oceans, Ottawa, Ontario K1A 0E6, (613) 990-0276 (Telephone), (613) 954-0807 (Facsimile).

#### Compliance and Enforcement

No new compliance mechanisms are being introduced by these amendments.

#### Contact

Ms. Iola M. Price, Director, Aquaculture and Oceans Science Branch, Department of Fisheries and Oceans, Ottawa, Ontario K1A 0E6, (613) 990-0275 (Telephone), (613) 954-0807 (Facsimile).

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## PROPOSED REGULATORY TEXT

Notice is hereby given that the Governor in Council, pursuant to section 43<sup>a</sup> of the *Fisheries Act*, proposes to make the annexed *Regulations Amending the Fish Health Protection Regulations*.

The proposed effective date of the Regulations is the date of registration thereof by the Clerk of the Privy Council.

Any interested person may make representations concerning the proposed Regulations within 30 days after the date of publication of this notice. All such representations must be addressed to Iola M. Price, Director, Aquaculture and Oceans Science Branch, Department of Fisheries and Oceans, 200 Kent Street, Station 1256, Ottawa, Ontario K1A 0E6, and cite the *Canada Gazette*, Part I, and the date of this notice.

December 19, 1996

MICHEL GARNEAU  
Assistant Clerk of the Privy Council

## REGULATIONS AMENDING THE FISH HEALTH PROTECTION REGULATIONS

## AMENDMENTS

1. Section 4 of the English version of the *Fish Health Protection Regulations*<sup>1</sup> is replaced by the following:

4. Subject to section 5, a local fish health officer for a province may issue, to a person who applies for one, an import permit that authorizes the person to import cultured fish or the eggs of wild fish into that province.

2. Sections 5 and 6 of the Regulations are replaced by the following:

5. No import permit shall be issued unless the person who applies for the permit has obtained a certificate and

(a) the certificate indicates that no disease or disease agent listed in Schedules II to IV was detected; or

(b) the local fish health officer is satisfied that none of the detected diseases or disease agents indicated on the certificate will be harmful to the conservation and protection of fish in the province of importation.

*Certificates*

6. A certificate required pursuant to section 5 is issued by a fish health official and

(a) certifies that the source of the fish was inspected in the approved manner; and

(b) indicates which, if any, of the diseases or disease agents listed in Schedules II to IV were detected during the inspection or inspections, as the case may be.

3. Item 6 of Schedule II to the Regulations is replaced by the following:

6. Enteric Redmouth Disease (*Yersinia ruckeri*)

## COMING INTO FORCE

4. These Regulations come into force on the date on which they are registered.

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[1-1-a]

<sup>a</sup> L.C. (1991), ch. 1, art. 12  
<sup>1</sup> C.R.C., ch. 812

## **NOTES ON THE PROPOSED INTERIM AMENDMENT TO FHPR**

### **Purpose of the amendment**

The purpose of this amendment is to:

1. Allow Local Fish Health Officers to approve transfers of eggs and fish between sites, even when disease agents of concern are detected at the source.
2. Allow Local Fish Health Officers to approve the transfer of disinfected eggs from source facilities or wild broodstocks that have only been tested for viruses.

It is not the intention of this amendment to supercede or cancel regional/provincial fish health protection policies.

## **GUIDELINES FOR IMPLEMENTING THE PROPOSED AMENDMENT TO THE FHPR**

### **A. IMPORTATION OF EGGS AND FISH**

#### **1. Importation of Disinfected Eggs of Wild and Cultured Fish**

Sources of eggs must be inspected by a Fish Health Official, and a Fish Health Certificate issued noting the presence or absence of filterable replicating viral agents, including their strains/serotypes if required by the recipient province. Viral agents include but are not limited to:

Viral Hemorrhagic Septicemia Virus (VHSV)  
Infectious Hematopoietic Necrosis Virus (IHNV)  
Infectious Pancreatic Necrosis Virus (IPNV).

The Local Fish Health Officer may, after reviewing the Fish Health Certificate, issue an Import Permit if the import of such eggs will not result in the introduction of a viral agent or strain/serotype of a viral agent listed above and not already known to occur in the receiving province.

In the absence of fish health data indicating the presence of selected pathogens, Local Fish Health Officers may designate areas in a province as free of selected pathogens, even though the pathogens have been detected in other parts of the province. Local Fish Health Officers may reject applications to import eggs, from a source that tested positive for a selected pathogen, into an area designated free of the specific pathogen.

The eggs must be accompanied by an Import Permit issued by the Local Fish Health Officer in the receiving province.

Sources of disinfected eggs need only provide information on virus testing when applying for Import Permits. If eggs are not disinfected, information on all diseases or disease agents listed in Schedule II of the Regulations (i.e. including bacteria and parasites) must be provided. Eggs must be disinfected at the source and in the receiving facility.

## 2. Importation of Live Cultured Fish

Sources of live cultured fish must be inspected by a Fish Health Official, and a Fish Health Certificate issued noting the presence or absence of diseases/disease agents listed in Schedule II of the FHPR, including:

Filterable replicating viral agents, including their strains/serotypes if required by the receiving province. Viral agents include but are not limited to :

Viral Hemorrhagic Septicemia Virus (VHSV)  
Infectious Hematopoietic Necrosis Virus (IHNV)  
Infectious Pancreatic Necrosis Virus (IPNV).

*Aeromonas salmonicida*  
*Yersinia ruckeri*

*Myxobolus cerebralis*  
*Ceratomyxa shasta*.

The Local Fish Health Officer may, after reviewing the Fish Health Certificate, issue an Import Permit if the import of such fish will not result in the introduction of a disease agent or strain/serotype of a disease agent listed above and not already known to occur in the receiving province.

In the absence of fish health data indicating the presence of selected pathogens, Local Fish Health Officers may designate areas in a province as free of selected pathogens, even though the pathogens have been detected in other parts of the province. Local Fish Health Officers may reject applications to import live cultured fish, from a source that tested positive for a selected pathogen, into an area designated free of the specific pathogen.

If fish at the source are clinically ill with Schedule II diseases or disease agents , no fish should be transferred until the disease episode is controlled.

The fish must be accompanied by an Import Permit issued by the Local Fish Health Officer in the receiving province.

### 3. Dead, Uneviscerated Cultured Fish

Sources of dead, uneviscerated fish must be inspected by a Fish Health Official, and a Fish Health Certificate issued noting the presence or absence of disease agents listed in Schedule III, including :

Viral Hemorrhagic Septicemia Virus (VHSV)  
*Myxobolus cerebralis*

The Local Fish Health Officer may, after reviewing the Fish Health Certificate, issue an Import Permit, if the import of such dead uneviscerated cultured fish will not result in the introduction of a disease agent or strain/serotype of a disease agent listed above and not already known to occur in the receiving province.

In the absence of fish health data indicating the presence of selected pathogens, Local Fish Health Officers may designate areas in a province as free of selected pathogens, even though the pathogens have been detected in other parts of the province. Local Fish Health Officers may reject applications to import dead uneviscerated cultured fish, from a source that tested positive for a selected pathogen, into an area designated free of the specific pathogen.

## **B. REQUIREMENTS TO OBTAIN A FISH HEALTH CERTIFICATE**

### 1 Aquaculture Facilities

An existing facility, with eggs and/or fish of unknown pathogen status, must have four inspections at approximately six-monthly intervals over a period of not less than 18 months before a Fish Health Certificate may be issued. Aquaculture facilities wishing to export disinfected eggs only, need only be tested for viral agents.

If eggs or fish are transferred from another source within a province, the Fish Health Certificate for the receiving facility will be changed to reflect the health status of the source (if disinfected eggs are transferred, only the status of Schedule II viruses will be changed for the receiving facility). If the health status of the source is unknown, that facility will have to re-commence four inspections at approximately six-monthly intervals to obtain a new Fish Health Certificate.

It will be the responsibility of the owner of the receiving facility to notify the Fish Health Official about any change in status of his/her Fish Health Certificate. The Fish Health Official will provide an amended Fish Health Certificate, copied to the National Registry of Fish Diseases. The expiry date on the amended Certificate will be the same as that on the voided Certificate.

A new production facility located on an isolated water supply free of all species of fish and starting with stocks from a source with a valid Fish Health Certificate, can obtain a Fish Health Certificate after only one inspection. The Fish Health Certificate must reflect the disease agent profile of the source facility, plus the results of the one inspection at the receiving facility.

Once a facility has received a Fish Health Certificate, two annual inspections at approximately six-monthly intervals are required to maintain the Fish Health Certificate.

2 Eggs of Wild fish

Wild adult broodstock must have a record of two consecutive annual inspections before a Fish Health Certificate may be issued. Where it is not possible to sample the same populations in consecutive years (e.g Pacific salmon), broodstock from the same area of the river must be tested in the second year.

3 Upgrading Fish Health Certificates

Persons wishing to upgrade their Fish Health Certificate, from positive to negative for a specific pathogen, will be required to:

- Implement a program to eliminate the specific pathogen(s).
- Have four consecutive negative [for the specific pathogen(s)] inspections over a minimum period of 18 months.

**C. ADDITIONAL DISEASE TESTING**

A Local Fish Health Officer in a receiving province may request additional testing in accordance with regional/provincial policies or regulations to determine the strain and/or nucleic acid profile for a pathogen detected at a source facility, if there is information indicating that a new strain or serotype of a pathogen listed in Schedules II, III and IV could be introduced to the receiving province with a shipment of eggs or fish.

If additional testing is required by the Local Fish Health Officer, then certain conditions of testing must be complied with. Additional testing must be coordinated and/or supervised by a Fish Health Official or Local Fish Health Officer. Tests conducted shall have a level of sensitivity and specificity comparable with other routine regulatory diagnostic tests, and test procedures must be accessible to both public and private laboratories doing health certification inspections. Additional testing must be conducted by a credible

diagnostician, in a timely fashion, and with due regard to provisions for continuity of evidence.

**D. REVIEW PROCEDURE (New)**

Where the Local Fish Health Officer does not issue an Import Permit under the Fish Health Protection Regulations, the Local Fish Health Officer must provide written reasons for rejecting an application.

The applicant may request a review of a rejected application by applying to the Assistant Deputy Minister (ADM), Science, Department of Fisheries and Oceans, for a review of the decision. The application should be in the form of a letter to the ADM addressed to the attention of the National Registrar. The application must be delivered to the National Registrar within 30 days of receipt of the decision by the Local Fish Health Officer not to issue the Import Permit.

The application must include a copy of the original application for the Import Permit, a copy of the reason(s) given for not issuing the permit, and the reason(s) for requesting a review. Any additional information which the applicant wishes to submit in support of the application must be included with the review request. The applicant must demonstrate that the decision not to issue the Import Permit was inconsistent with the Regulations.

The National Registrar will coordinate the review process, establish a Review Board, participate as an advisor on the Review Board, hold the documentation to be forwarded to the ADM until all of the documentation is assembled and forward the documentation to the ADM.

The National Registrar will forward a copy of the application to the appropriate DFO Regional Director General and responsible provincial agencies with a request that the Regional Director General and responsible provincial agencies submit written comments to the Review Board with respect to the application. They should respond within 21 calendar days of receiving the application from the National Registrar.

The Review Board will be established within 5 calendar days of receipt of the application and will consist of the National Registrar and three independent persons selected by the National Registrar in consultation with the applicant and the responsible federal and provincial agencies. The role of the National Registrar on the Review Board will be advisory. The National Registrar will not vote on the recommendation(s) set out in the report prepared by the Review Board.

The role of the Review Board in the review process is advisory. The Review Board provides a report, including its recommendation(s) to the ADM. The ADM makes the final decision whether an Import Permit should be issued.

The Review Board may request written submissions from anyone who, in the opinion of the Review Board, should be consulted.

The Review Board will prepare a written report to the ADM within 30 days of the Review Board being established. In exceptional circumstances, the Review Board may request the applicant for a reasonable extension of time to complete its report. The report will contain the recommendation(s) of the Review Board with respect to the application.

Subject to the Access to Information Act and the Privacy Act as amended from time to time, the Review Board will provide copies of the documentation to be forwarded to the ADM, to the applicant. The documentation will include a copy of the Review Board's report and recommendation(s) and the written submissions of the Regional Director General and the responsible provincial agencies and other written submissions obtained by the Review Board as part of its review.

The applicant will have 10 calendar days from receipt of the documentation to forward a written submission setting out comments. The submission should be in the form of a letter to the ADM addressed to the attention of the National Registrar.

The National Registrar will forward the documentation to the ADM for review and a final decision. The documentation will include the applicant's application for review (including any supporting documentation filed with the application), the written submissions of the Regional Director General and responsible provincial agencies and any other written submissions obtained by the Review Board in the course of its review, the Review Board's report and recommendation(s), and the applicant's written submission setting out comments.

The ADM will make the final decision on whether to issue an Import Permit and will provide the decision and the reason(s) for the decision in writing within 30 calendar days. The decision will be forwarded to the applicant, copied to the Regional Director General, responsible provincial agencies and the Local Fish Health Officer, through the National Registrar.

**Note:** If a review relates to a shipment of eggs, any juvenile fish that may hatch from these eggs before the completion of the review procedure will be subject to live fish health requirements of the amended Fish Health Protection Regulations and other related regulations and policies.

### SUMMARY OF REVIEW PROCEDURE

ACTION	ALLOTTED TIME PERIOD	EXPECTED ELAPSED TIME AFTER LFHO SENDS REJECTION LETTER TO APPLICANT
LFHO sends letter to applicant rejecting application to import eggs or fish	-	Day zero
Applicant sends review to National Registrar	30 days	Day 30
National Registrar establishes Review Board	5 days	Day 35
DFO Regional Director General and responsible provincial agencies submit comments on application to National Registrar	21 days	
Review Board prepares report and recommendations for ADM	30 days	Day 65
Applicant provides comments to National Registrar on Review Board's report and recommendations	10 days	Day 75
ADM provides written decision on appeal to applicant	30 days	Day 105

Aquaculture and Ocean Science Branch  
Department of Fisheries and Oceans  
Ottawa

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